

# Cleaning and Disinfecting

## Stainless Steel Cleaner & Polish

- Powerful, non-corrosive cleaner for all stainless steel surfaces
- Super concentrated biodegradable formula cleans and polishes in one simple step
- Removes fresh or dried blood and tissue
- Ideal for autopsy tables, grossing stations and all stainless steel instrumentation including:
  - Biomedical Instrumentation
  - Lighting Fixtures
  - Freezers and Coolers
  - Refrigerators
  - Body Trays and Carriers
- Fresh citrus fragrance
- 32 oz. (.95L) spray bottle



Product No. BE048

## Deodorx™

Developed especially for use by Medical Examiners, Deodorx™ eliminates offensive odors emitted from decomposed bodies within seconds.

From the discovery scene to the morgue, Deodorx can be used in every phase of your work. Deodorx not only removes odors from clothes, vehicles, body bags, refrigerators and other items, it keeps the odor from permeating the lab environment.

New, improved formula and completely water soluble, odorless, colorless, biodegradable, non-flammable and non-corrosive. Contains no CFC's, formaldehyde or hazardous biocides; will not attack metals, plastics, or fabrics and emits no toxic fumes. Sold in one gallon (3.8L) containers.

Product No. BE040



## Lysol® I.C. Disinfectant Concentrate

Impressive anti-microbial activity - Disinfects in ten minutes!  
Concentrated formula provides low cost disinfection in the pathology/necropsy areas. One gallon (3.79L) disinfects approximately 140,000 square feet of surface at a recommended 1:50 dilution in a hospital setting, 1:256 dilution, broad spectrum use. Destroys a wide range of pathogenic bacteria, viruses and fungi. Meets OSHA's Bloodborne Pathogens Standard for decontamination of hard, non-porous environmental surfaces such as stainless steel. EPA Registration No. 47371-129-675.



Product No. BE022

## Deodorx Plus™

Also available is Deodorx Plus - Super cleaner deodorizer and degreaser for the morgue, lab and pathology equipment, including grossing stations, refrigerator, body bags, cabinets, tables, sinks and more. Sold in one gallon (3.8L) containers.

Product No. BE042



## EZ-Zyme Instrument Cleaner



Simply open the packet of concentrate. Dilute into a gallon of warm tap water and place your hard to clean instrument into the solution. No hand scrubbing required. All tissue debris will dissolve leaving every crack and crevice of your dissection instruments looking brand new. This product is not a disinfectant, just one of the best enzyme cleaners on the market. Easy to use and, saves storage space. Box contains 32 packets.

Product No. BE100



## Spray Lube

Spray on all moving parts of instruments after cleaning and before autoclaving. Let instrument drain. Rinsing and wiping are unnecessary. Instruments must be autoclaved. Contains 8 fluid oz.

Product No. BE060

1

2

3

4

5

6

7

8

9

10

11