



Disinfectants and Cleaners

*Whatever the Application,
Mopec Has the Right
Formula for You!*

**Specially
Formulated
For:**

- ✓ **Pathology**
- ✓ **Morgue**
- ✓ **Histology**
- ✓ **Necropsy**



Kills TB and HIV-1



*Ideal for Pathology and
Morgue Surfaces*



*No More Scrubbing
Dissecting Instruments*



*Super Powerful – Clean
and Polish in One Step*



*A wipe that contains
no bleach and
is non-phenylphenol*



*Effortlessly clean and
polish all your lab equipment
at the same time*

Mopec The Better Product People

Call Us at: **1-800-362-8491** or Visit Us on the Web at **www.mopec.com**

PATH CLOUD

This powerful foaming action cleans, disinfects and deodorizes in one easy operation. EPA approved, Path Cloud has been proven effective against HIV-1, Herpes Simplex 2, Staphylococcus, Pseudomonas, Canine Parvo, Tuberculosis and other infectious diseases. Also destroys a wide range of pathogenic bacteria, viruses and fungi. The tuberculocidal meets OSHA's Bloodborne Pathogens Standard for decontamination of hard, non-porous environmental surfaces. Sold in a 32 oz. spray bottle. Ideal for heavily soiled pathology surfaces including:

- Autopsy and dissection tables
- Grossing stations
- Morgue/cutting room sinks
- Laboratory counters/backsplashes
- Dissection Boards
- Body trays

EPA # 5747-MI-1 EPA Registration No. 1839-83-7495

Order Number: BE045



BENCH SPRAY & WIPE

Specially formulated multi-purpose decontaminating cleaner safe for all work surfaces and equipment including: histology cutting stations, biomedical instrumentation, morgue/cutting room floors, computer work stations, laboratory telephones and dictation equipment. Ideal for bench tops and tables, refrigerators, freezers and laboratory appliances. Bench Spray & Wipe is EPA approved and has proven effective against HIV-1, Staphylococcus, Herpes Simplex 2, , Pseudomonas, Canine Parvo, Tuberculosis and more. This specially formulated tuberculocidal meets OSHA's Bloodborne Pathogens Standard for decontamination of hard, non-porous environmental surfaces. Mild lemon fresh fragrance, sold in a 32 oz. spray bottle.

EPA # 5747-MI-1 EPA Registration No. 1839-83-74952

Order Number: BE047



EAZY SOAK

Enzyme instrument presoak and cleaner. Our unique enzyme instrument cleaner works without scrubbing or pre-rinsing. Eazy Soak digests all tissue, blood and exudative materials. Perfect for use on hard-to-clean scissors, forceps, serrated cutting blades and scalpel handles. Clean your lab instruments in three easy steps. Simply add 2 oz. Eazy Soak to one gallon of hot tap water, add the soiled instruments and walk away. Your instruments will be clean without the hard work. Concentrated liquid sold in 32 oz. containers.

Order Number: BE046



STAINLESS STEEL CLEANER AND POLISH

Powerful, non-corrosive cleaner safe for all stainless steel surfaces. Super concentrated biodegradable formula cleans and polishes in one simple step. Removes fresh or dried blood and tissue easily. Ideal for autopsy tables, grossing stations and all stainless steel instrumentation including: Biomedical Instrumentation, Bone Saws, Cutting Surfaces, Stainless Steel Freezers, Refrigerators, morgue coolers, body trays and carriers. Fresh citrus fragrance, sold in 32 oz. spray bottle.

Order Number: BE048



PATHWIPES

One step cleaning and disinfecting of all pathology and necropsy and lab surfaces. EPA approved. Pull-out disposable towelettes. Effective against HIV-1, Herpes Simplex 2, Pseudomonas, MRSA, Staphylococcus, Tuberculosis and Canine Parvo. Ideal for lab telephones, keyboards, writing instruments and splash shields. Perfect for cleaning and disinfecting laboratory equipment, lab refrigerators, instruments, doorknobs, computer mouse, rulers and microscopes. Don't just clean those last minute blood spots and smears - disinfectant at the same time. 160 pre-measured wipes, 6" x 6-3/4" in each canister.

EPA # 5747-MI-01 EPA Registration No. 1839-174-74952

Order Number: BE036



STAINLESS STEEL CLEANING AND POLISHING WIPES

Brings out the natural luster in your equipment while leaving a protective covering that repels water, fingerprints and dust. Especially designed to be used on all lab equipment such as workstations, cabinets, refrigerators and storage areas. Stainless Steel Cleaning and Polishing Wipes can restore a brilliant "like new" shine to your entire lab. Removes stains, smudges and blood while leaving behind a highly polished, shiny finish. Clean-up is quick and thorough. Simply pull wipe from canister and with moderate pressure clean soiled area. 105 pre-measured wipes, 7" x 10" in each canister. Allow to dry for 5 minutes.

Order Number: BE039



For a complete catalog detailing our products and services, please call or visit our web site



The Better Product People

21750 Coolidge Highway • Oak Park, MI 48237
Toll Free: 800-362-8491 • Tel: 248-291-2040 • Fax: 248-291-2050

www.mopec.com